


Faster Triage, Prompt Actions
Because Time Matters.



Automated Calculation of Cerebral Hemodynamics Heuron CTP

- ✓ CT Perfusion Analysis
- ✓ Faster Decision-Making
- ✓ Customized Viewer Options



Analyzes cerebral hemodynamics for a fast selection
of patients in need of emergency treatment



Workflow Improvement

Automatically analyzes CT perfusion and provides clear, informative hemodynamics with color-coded maps.
(CBF, CBV, MTT, TTP, Tmax)



Precise Information

Provides essential information for treatment decisions, including Tmax Volume (>6sec, Penumbra) and rCBF Volume (<30%, Infarct Core), and Mismatch Map.



AIF/VOF Automation

Automatically generates AIF/VOF graphs and allows user customization for reanalysis, enabling more precise evaluation of cerebral blood flow and tissue characteristics.



台北市信義區忠孝東路五段550號14樓

T. +886-2-2346-3468 | E. kcpigi09458008@gmail.com | H. www.iheuron.com

© 2025. Heuron Co., Ltd. All rights reserved.

Heuron-PB-CTP-02E (03)



Scan to Connect and
Discover More